Three-dimensional gantries EXCL





Key features

At a glance

The three-dimensional gantry EXCL is ideal, for example, for analytical processes where sample vials are to be opened and liquid samples pipetted in the same three-dimensional gantry. The optional second Z-axis enables two front units to move independently of each other, such as the rotary gripper module EHMD and the smart pipette DHOP.

Application examples with further Festo components

- Automated sample preparation for the analysis of sample vials (e.g. clinical analysis of patient samples or material analysis of samples for quality assurance)
- Picking up sample vessels from a rack and rotating them

 using the rotary gripper module EHMD
- Taking a liquid sample and transferring it to a microwell plate for further analysis
- using the smart pipette DHOP or pipette head DHOE
- Dispensing reagents into the microwell plate – using the dispense head VTOE
- Implementing rinsing procedures in the ELISA process
- using the dispense head VTOI

Operating mode



EXCL is a three-dimensional gantry for installation in devices such as analytical table top units and is optimised for life sciences applications. The working area (X,Y) can be configured up to 1000 mm x 700 mm in 1-mm steps. The stroke of each of the Z-axes can be selected separately from 50, 100, 150 and 200 mm. Configuration as a planar surface gantry X,Y without Z-axes is also possible.

Key features



Key features



001	Series	008	Stroke of the Z-axis 2	
EXCL	Three-dimensional gantry		No Z-axis	
		50	50 mm	
002	Size	100	100 mm	
15	15	150	150 mm	
	· · ·	200	200 mm	
003	Stroke of the X-axis [mm]			
200	200	009	Guide	
400	400	GF	Plain bearing	
1000	1000	010	Motor type	
004	Stroke of the Y-axis [mm]	ST	Stepper motor ST	
200	200			
300	300	011	Motor attachment position	
700	700	L	Left	
1		R	Right	
005	Number of Z-axes			
	No Z-axis	012	Cable length [m]	
Z1	One Z-axis	0.2	0.2	
Z2	Two Z-axes	1	1	
006	Version of the Y-axis	2	2	
	Standard	013	Energy chain	
YS1	1 slide	EK	With	
YS2	2 slides			
	<u>_</u>	014	Controller	
007	Stroke of the Z-axis 1		None	
	No Z-axis	M6	For six motors	
50	50 mm			
100	100 mm	015	Position sensing	
150	150 mm	AN	With integrated zero position sensing	
200	200 mm			·•

Peripherals overview



Acces	Accessories							
Туре		Description	→ Page/Internet					
[1]	Three-dimensional gantry EXCL	Electric	7					
[2]	Servo drive CMXL	PCB with 6-axis motion controller for direction integration into the device	15					
[3]	Adapter kit EHAM	 For mounting the gantry via the end cap (shown on one side by way of example) 4 included in the scope of delivery 	15					

- 📲 - Note

The stepper motors, motor cables and connecting cables to the reference switches are included in the scope of delivery of the three-dimensional gantry. The servo drive is optional. i

Datasheet



General technical data		
Design		Three-dimensional gantry system
Guide		Plain-bearing guide
Stroke of the		
X-axis	[mm]	200 1000
Y-axis	[mm]	200 700
Z-axis	[mm]	50, 100, 150, 200
Min. stroke limit	[mm]	50
Max. stroke limit	[mm]	1000
Rated load ¹⁾	[kg]	1.5 (with one Z-axis)
	[kg]	2 (in total, with 2 Z-axes)
Max. process force ²⁾	[N]	15
Motor holding torque	[Nm]	0.5 1.4
Rated voltage		
X-axis	[V]	1.68
Y-axis	[V]	3.2
Z-axis	[V]	3.2
Rated current		
X-axis	[A]	4.2
Y-axis	[A]	2
Z-axis	[A]	2
Max. acceleration		
X-axis	[m/s ²]	1
Y-axis	[m/s ²]	1
Z-axis	[m/s ²]	0.2
Max. speed		
X-axis	[m/s]	0.4
Y-axis ³⁾	[m/s]	0.5
Z-axis	[m/s]	0.07
Max. search speed velocity		
X-axis	[m/s]	0.03
Y-axis	[m/s]	0.03
Z-axis	[m/s]	0.015
Max. creep speed		
X-axis	[m/s]	0.003
Y-axis	[m/s]	0.003
Z-axis	[m/s]	0.003
Repetition accuracy	[mm]	±0.1
Mounting position		Horizontal
Type of mounting		With female thread M3

1) Rated load = tool load (attachment components) + payload

2) Perpendicular to the working plane, at standstill

3) Tested at: v=0.4m/s and two Z-axes: Additional weight: Z1=0.5kg and Z2=1.5kg

Datasheet

Operating and environmental conditions

Operating and environmental conditions		
Degree of protection		IP10
Ambient temperature ¹⁾	[°C]	+10+40
Storage temperature	[°C]	-10+60
Relative humidity	[%]	0 90 (non-condensing)
Noise level	[dB(A)]	79.2
Duty cycle	[%]	100
Vibration resistant		Transport inspection to FN 942017-4 and EN 60068-2-6
Shock resistance		Shock test with severity level 1 to FN 942017-5 and EN 60068-2-27
КС		KC EMC
CE marking (see declaration of conformity) ²⁾		To EU EMC Directive
UKCA marking (see declaration of conformity)	2)	To UK EMC regulations

1) Note operating range of proximity sensors and motors

More information www.festo.com/catalogue/excl → Support/Downloads

Materials		
X/Y/Z profile		Aluminium
Housing		Aluminium
Guide		РОМ
X-axis drive pinion		Aluminium
Y-axis drive pinion		РОМ
Clamping pinion		Aluminium
Toothed belt		Polychloroprene with glass cord
Gear rack		Aluminium
LABS (PWIS) conformity		VDMA24364 zone III
Note on materials		RoHS-compliant
Toothed belt for X-axis		
Pitch	[mm]	2
Width	[mm]	8
Feed constant		
X-axis	[mm/rev]	64
Y-axis	[mm/rev]	37.7
Z-axis	[mm/rev]	5
Weight [kg]		
Product weight (without rated load a	nd controllers)	
Stroke: x = 200 mm; Y = 300 mn	n; 1x z = 50 mm	5.45
Stroke: x = 400 mm; Y = 700 mn		9.12
Base weight		
With 1 slide		3.68
With 2 slides		4.66
Additional weight per 50 mm stroke		
X-axis		0.14
Y-axis with 1 slide		0.09
Y-axis with 2 slides		0.13
Weight		
Adapter kit EHAM		0.35

Dimensions

Dimensions

EXCL-15-...-Z1-...-L





+ = plus stroke length

Туре	B1	B2	B3	B4	B5	B6	H1	H2	H3	H4
[mm]										
EXCL-15	202	98.3	147	77	200 700	24	214	144	85.3	50
Type [mm]	H5	H6	L1	L2		3	L4	L5	L6	L7
EXCL-15	14	50 200	219	215.	5 16	5.5	200 1000	50	25	5

Download CAD data → <u>www.festo.com</u>

Dimensions

Dimensions





Download CAD data → <u>www.festo.com</u>

Dimensions



EXCL-15

43

214

214

Dimensions

Dimensions

۲9 ا

TЗ



Туре

[mm]

Туре

[mm] EXCL-15 H7

29

H8

25

H9

20.5

H10

12

H11

4

L8

20

L9

5.5

T1

6

T2

1.3

T3

1.3

T4

6.5

T5

1.3

T6

2.6

T7

1.6

EXCL-15

Dimensions

Dimensions

EXCL-15-...-YS... (interface)





Туре	B12	D11	D12	L10	T8	Т9
		Ø				
[mm]	±0.05	H7		±0.05	max.	
EXCL-15	40	5	М3	40	6	1.3

Download CAD data \rightarrow <u>www.festo.com</u>

Ordering data – Modular product system

Ordering table					
Size		15	Conditions	Code	Enter cod
Module no.		8159775			
Product type		EXCL series L		EXCL	EXCL
Size		15		-15	15
Stroke of the X-axis	[mm]	200 1000			
Stroke of the Y-axis	[mm]	200 700			
Number of Z-axes		No Z-axis			
		1 Z-axis		-Z1	
		2 Z-axes		-Z2	
Y-axis version		Standard			
		1 slide	[1]	YS1	
		2 slides	[1]	YS2	
Stroke of the Z-axis 1		No Z-axis			
	[mm]	50	[2]	-50	
	[mm]	100	[2]	-100	
	[mm]	150	[2]	-150	
	[mm]	200	[2]	-200	
Stroke of the Z-axis 2		No Z-axis			
	[mm]	50	[3]	-50	
	[mm]	100	[3]	-100	
	[mm]	150	[3]	-150	
	[mm]	200	[3]	-200	
Guide		Plain-bearing guide		-GF	-GF
Motor type		Stepper motor ST		-ST	-ST
Motor attachment position		Left		۰L	
		Right		-R	
Cable length ¹⁾		0.2 2 m			
Energy chain		with		-EK	-EK
Controller		None			
		For 6 motors		-M6	
Position sensing		With integrated zero position sensing		-AN	

[1] YS1, YS2 Mandatory with number of Z-axes "no axis"; Not with number of Z-axes "1 axis"; Not with number of Z-axes "2 axes"

[2] **50, 100, 150, 200** Mandatory with number of Z-axes "1 axis"; Mandatory with number of Z-axes "2 axes"

[3] 50, 100, 150, 200 Mandatory with number of Z-axes "2 axes"

Selection of cable length

1) The cable length can be selected using the modular product system. This specification relates to the output of the energy chain at the X-axis (dimension L) and describes the minimum length by which the cables protrude.



Ordering data – Three-dimensional gantry

Ordering data – Three-dimensional gantry										
	Stroke of the	Stroke of the	Number of	Stroke of the	Part no.	Туре				
	X-axis	Y-axis	Z-axes	Z-axes						
	200	300	1	100	8184925	EXCL-15-200-300-Z1-100-GF-ST-R-1-EK-M6-AN				
	400	700	2	200	8184926	EXCL-15-400-700-Z2-200-200-GF-ST-R-1-EK-M6-AN				

Three-dimensional gantries EXCL

NEW

Accessories

Adapter kit EHAM

Material:

Plates: Anodised aluminium Screws: Galvanised steel Centring sleeves: Stainless steel RoHS-compliant





Β1

[1] 2x ISO 14579-M3-16

Dimensions an	Dimensions and ordering data											
For type	B1		B2	E	33	B4		D1	H1		H2	
[mm]								Ø				
EXCL-15	40		21		7	26		6	14		6	
			1			1					_	
For type	L1	L2	L3	L4	L5	L6	L7	Weight	Part no.	Туре		
[mm]								[g]				
EXCL-15	100	75	50	30	25	72	14	355	8185812	FUAAA	MA-E23-15	

Servo drive

Designation	Description	Nominal voltage Load supply [V DC]	Number of digital logic inputs/outputs	Part no.	Туре
	 PCB with 6-axis motion controller for the direct control of the EXCL's stepper motors. Pre-configured for the three-dimensional gantry EXCL (motor and axis data/gantry size, homing) 6 axes, 4.5 A per motor, open-loop, with inputs for the EXCL microswitch for homing Programming with G-code For direct integration in devices, not for H-rail mounting 	24	9/19	8168593	CMXL-ST6-NE-C4-1-MP